



Stability Statement

Part Number: [33814-Glu](#)

ProductDescription: Glucose Test Strips 0-3000mg/dL

The test strips should be stored in the closed bag in a cool, dry place out of direct sunlight. The ideal condition is 25C (77F). A range of 20C (68F) - 30C (86F) is acceptable. Test strips should not be stored in a refrigerator or freezer. Test strips should not be stored where they can be exposed to heat (for example, inside a vehicle or shed), damp or high humidity environments. The bag these strips are shipped in is amber to prolong their shelf life and reactivity.

When stored properly, the test papers have a minimum shelf-life of 24 months from the date of manufacture. Under proper storage conditions the strips will continue to perform properly beyond the expiration date.

Indigo Instruments

169 Lexington Crt. Unit i

Waterloo, ON N2J 4R9

Phone: (519) 746 4761

Email: info@indigoinstruments.com

Web: www.indigoinstruments.com